

EN

IDEAL

MADE IN ITALY

Extractive Surgery



ASA DENTAL
MAKE PEOPLE SMILE

IDEAL by Asa Dental

GENTLE & PRECISE



ATRAUMATIC EXTRACTION

Remove the tooth with confidence through an atraumatic luxation



ERGONOMIC DESIGN

Take advantage of the handle design to apply only the necessary force and have total control



SHAPE AND SIZE

Operate with the most suitable tool for the case at hand, choosing between tips of different shapes and sizes

Dental care is constantly evolving to introduce those techniques and tools that ensure greater customer satisfaction. In case of invasive interventions, such as the removal of a tooth, outcome and patient satisfaction are closely linked to minimizing trauma. This is exactly the task that Asa Dental's new IDEALS perform precisely and effectively.

The IDEALS are atraumatic luxating instruments which, thanks to the sharp tip, gently and precisely cut the periodontal ligaments, saving unnecessary soft tissue injury and preventing fractures in the operation of extracting the tooth from the alveolus.

IDEAL Periotomes

With their sharp tips, these periotomes are designed to gently cut periodontal ligaments, preserving the surrounding gum and the bone.

To minimize the risk of fracture, the tooth is dislocated only after a careful expansion of the alveolus using the precise movements that the ergonomic handle allows.

0280-3R

Luxating Elevator,
Straight #3, Silver, 3 mm

0280-2R

Luxating Elevator,
Straight #2, Emerald, 2 mm

IDEAL Periotomes



0280-3C
Luxating Elevator,
Curve #3, Jet Black, 3 mm

0280-5R
Luxating Elevator,
Straight #5, Golden, 5 mm



IDEAL Periostomes

S0280-4

IDEAL Selection is more than a simple set of instruments. Cherry-picked by Asa Dental, the 4 included periostomes help you confidently overcome any obstacle typically associated with surgical extractions. The design and manufacturing of each instrument make luxation atraumatic to minimize treatment trauma and bleeding.

0280-5C

Luxating Elevator,
Curve #5, Jet Black, 5 mm

IDEAL Dual Edge Periotomes

These periotomes have a dual edge tip that facilitates penetration and reduces the risk of the instrument slipping. For this reason, they are particularly suitable in case of difficult extractions or broken teeth.

0280-5DF

Dual Edge Luxating Elevator,
#5, Gold, 5 mm

0280-3DF

Dual Edge Luxating Elevator,
#3, Silver, 3 mm



ORDERING INFO

IDEAL Periotes

code	0280-2R	Straight	#2	Emerald	2 mm
	0280-3R	Straight	#3	Silver	3 mm
	0280-3C	Curve	#3	Jet Black	3 mm
	0280-5R	Straight	#5	Golden	5 mm
	0280-5C	Curve	#5	Jet Black	5 mm

IDEAL Dual Edge Periotes

code	0280-3DF	Dual Edge	#3	Silver	3 mm
	0280-5DF	Dual Edge	#5	Golden	3 mm

IDEAL Selection

code	S0280-4				
includes:				0280-3R, 0280-3C, 0280-5R, 0280-5C	

HEAD OFFICE

Asa Dental S.p.A.
Via Valenzana di Sopra, 60
Bozzano
55054 Massarosa (Lu)
ITALY

Tel. +39 0584 938305
+39 0584 93363
+39 0584 938306
Fax +39 0584 937167

asadental@asadental.it
www.asadental.com

LOGISTICS CENTER

Via Sarzanese Valdera, 1331
Bozzano
55054 Massarosa (Lu)
ITALY

PRODUCTION PLANTS

Via delle Piagge, 1/A - Martia
55014 Capannori (Lu)
ITALY

Via Fierla, 13 - Z.I.
33085 Maniago (Pn)
ITALY

BRANCH OFFICES

Asa Dental USA Corp.
105 Maxess Road, Suite 124
Melville, New York 11747 - USA

Asa Dental CHINA CO. LTD.
Room 315, Huaying Building,
Zhong Xin Avenue
Tianjin Airport Economic Area
Tianjin - CHINA

Asa Dental S.p.A. (Branch)
Dubai South Business Center
Building C, 3rd Floor
Dubai Logistics City
PO Box 390667
Dubai - UAE

